

	A	B	C	D	E	F
1	Cat Num	Description	QTY	Price	Shape	Notes
2		<b>Leica accessories newer style</b>				
3	565-015	Leica LU 4/25 Epi Fluor Illum 4 cube	1	\$800	4	Modified for OptiGrid great shape
4	565-015	Leica LU 4/25 Epi Fluor Illum 4 cube	1	\$700	4	Modified for OptiGrid great shape
5						
6		<b>LEICA STEREO MICROSCOPE</b>				
7	445-617	Leica microscope carrier for focus drive	2	150	4	For MS5/MZ6,8,12
8	418-587	Leica coarse focus mount for TV zoom	1	200	4	works on 424-300
9	424-300	Leica Plain stand w B&W stage plate	1	200	4	works on 418-587
10						
11		<b>Leitz Wild Objectives 160/170mm TL</b>				
12	175-096	Wild H100 na 1.25 oel iris	1	293	4	not checked
13	519-249	Leitz 100x na 1.30 oel 170/.17	1	310	4	not checked lens looks good
14	519-249	Leitz 100x na 1.30 oel 170/.17	1	25	1	front lens broken for parts
15	519-382	Leitz Adapter negative lens Flu	1	100	4	not checked
16	519-419	Leitz 40x na .65 170/.17	1	210	4	quick check works
17	519-519	Leitz 25x Apochromat na .65 170/.17	1	330	4	quick check works .65WD
18	519-538	Leitz L32x na .40 phaco 1 170mm	1	275	4	quick check works 6.45 WD
19	519-646	Leitz 25x ach na .75 oel 170/.17flz	1	175	4	quick check works .36 WD
20	519-803?	Leitz ? Oel 100/fl oel d 1.20 obj insert	1	5	NOS	
21	559-002	Leitz UM 5x na .10 wd 1.4 aper	1	150	4	
22		<b>Leitz Eyepieces for 160/170mm TL</b>				
23	519-137	Periplan GF 10x EP fixed	4	100ea	4	
24		23.2 mm Dia for 170mm TL				18mm FOV
25	519-463	Periplan GF 12.5x M EP focus	3	100ea	4	18mm FOV
26		23.2mm for 170mm TL				
27	519-030	6x B EP fixed 18 mm FOV	2	100ea	4	
28		23.2mm for 170mm TL				
29	511-028	G 12.5x Large Focus 31.5 mm	2	150ea	4	<b>For FS stereo</b>
30	511-095	A 15x 31.5mm fixed	2	150ea	4	For FG1 FG2 Stereo
31		Periplan 10x Mess focus	1	100	3	photo reticle slight scratch
32	519-613	Periplan 10x small newer style	1	125	3	missing screw on eyeguard
33		23.2 mm dia for 160 mm TL				

	A	B	C	D	E	F
1	Cat Num	Description	QTY	Price	Shape	Notes
34	519-414	Periplan 10x M Focusable	4	100ea	4	15 mm FOV
35		23.2 mm dia for 170mm TL				
36	519-170	Periplan 10x fixed 14mm FOV	1	100	4	w reticle
37		23.2mm dia for 160mm TL				
38	519-613	Periplan 10x o EP fixed	1	100	4	
39		23.2mm dia for 160mm TL				
40						
41		<b>Wild Eyepieces and stereo aux lens</b>				
42	346-922	15x/17 Eyepiece focus for M400	1	150	4	w photo reticle
43	360-208	Wild 15x/17 Eyepiece fixed for M400	1	150	4	For M400 pair w 346-922
44	198-130	10x K EP	2	75ea	4	
45	175-095	6x K EP	1	75	4	Slight dinge
46	131-515	10x PK EP	1	75	4	
47	<b>367-898</b>	<b>Wild .4x attachable f M400</b>	1	380	NOS	
48	<b>383-978</b>	<b>Wild .5x attachable WD 188 f M400</b>	1	380	NOS	
49						
50		<b>Leitz Parts and accessories</b>				
51	XXXX-1	Leitz Phase annular diaphragm	1	125	NOS	for 513-140 Phaco No. 402a cond Pos 1
52	XXXX-2	Leitz Phase annular diaphragm	1	125	NOS	for 513-140 Phaco No. 402a cond Pos 2
53	XXXX-3	Leitz Phase annular diaphragm	1	125	NOS	for 513-140 Phaco No. 402a cond Pos 3
54	513-547	Leitz Standard cond base LK w aperat	1	200	NOS	For Laborlux 12 type
55	512-420	Leitz .90 Ach condenser top	2	120ea	3	screw into condenser 512-081
56	513-450?	Leitz Side by side dual view	1	480	3	dark grey
57	514-073	Leitz Heat sink	1	93	NOS	
58	514-602	Leitz Illum House 6v 20w No 20	1	289	NOS	in original box for Lab 12 excludes socket
59	514-674	Leitz 6v 20w socket for LH 20	4	190 ea	NOS	for 514602 illumin house 20 for Lab 12
60	956-202	Leitz HM Lux metal back	1	45	NOS	For epi illuminator plate
61		Leitz Screws for bridge	2	1.00 ea	NOS	
62	051-406	Leitz Feet for LM Lux	3	8ea	NOS	
63		Leitz Small brass bushing for 1512 microto	1	25	NOS	
64	514-027	Leitz Heat filter for LH 250 53mm 2mm	1	20	NOS	
65		Leitz Slide holder	1	75	NOS	curved finger part

	A	B	C	D	E	F
1	Cat Num	Description	QTY	Price	Shape	Notes
66		Leitz HM Lux metal back	1	10	NOS	Solid metal back plate
67		Leitz Condenser dovetail for Ortholux II	1	170	NOS	
68		Leitz HM Lux metal field diaphragm	2	20 ea	NOS	
69		Leitz stage screws	4	3 ea	NOS	
70	513-467	Leitz Light ring 1 S1.1 condenser insert	1	121	NOS	phase ring insert for 513-467 univ cond
71	020-422.521-011	Leitz slide holder finger clip spring	1	125	NOS	
72	026-393-811-015	Leitz slide holder finger clip spring	1	125		
73	956-202	Leitz metal back plate	1	175	4	Light color
74		Leitz Illum mount plate	1	150	4	light color for SM Lux
75		Leitz Nosepiece clips	2	120ea	NOS	
76		Leitz Blue filter plastic	2	15ea	4	For Labolux 11
77		Leitz Green filter glass	1	20	4	For Labolux 11
78	513-362	Leitz oel 100/f/ 95 + D 1.20 funnel stop	1	55	4	for achromat obj 100/1.30
79		<b>Wild Repair parts</b>				
80	23485	Wild Phase ring 31 mm	2	120ea	NOS	
81	162-555	Wild Screw	5	16 ea	NOS	
82						
83	026-360-001-058	Wild Spring	1	17	NOS	
84	1042-81	Wild thumb wheel	2	21 ea	NOS	
85	108-122	Wild Stereo stage insert ground glass	1	25	4	80mm dia 50mm f M3/5/7
86	125-250	Wild prism	1	75	NOS	
87	125-276	Wild Puffer	5	14 ea	NOS	
88	125-277	Wild Screw	5	18 ea	NOS	
89	<b>125-284</b>	<b>Wild Pinion</b>	<b>1</b>	<b>110 ea</b>	<b>NOS</b>	<b>for stand serial # 97545</b>
90	126-126	Wild Daylight filter 32mm 2mm	2	12ea	NOS	
91	127-735	Wild Clamping band	1	35	NOS	
92	150-ps11	Wild Switch	1	9	NOS	
93	153-308	Wild Nosepiece clip for M-40	4	35ea	NOS	
94	162-555	Wild Screws	5	12ea	NOS	
95	167-720	Wild Plugs	3	10ea	NOS	
96	170-962	Wild bolts that hold rack on	4	10ea	NOS	f stand serial #below 94335
97	175-278	Wild Spring	3	12 ea	NOS	

	A	B	C	D	E	F
1	Cat Num	Description	QTY	Price	Shape	Notes
98	177-577	Wild R knob threaded	6	75ea	NOS	f stand serial #below 94335
99	178-859	Wild Focusable condenser holder	1	135	NOS	
100	108-122	Wild Stereo stage insert acrylic diffused	1	15	3	80mm dia 50mm thick f M3/5/7
101	184-506	Wild Eyepiece Guard hard plastic	2		NOS	Fit on wild 15x ep that I have
102	184-657	Wild Illumination lens	3	110ea	NOS	1 lens
103	184-659	Wild Illumination lens	1	95	NOS	
104	185-242	Wild screw	2	4 ea	NOS	
105	186-838	Wild 50x50 1/4 wave sheets	1	45	NOS	
106	194-674	Wild Screw	6	11ea	NOS	
107	195-181	Wild MPS 0.8x lens	1	65		
108	195-196	Wild MPS 1.0x lens	1	55	NOS	
109	196-658	Wild Inverted Stage insert	1	45	NOS	standard size slide holder
110	212-428	Wild cap and Screw 1(ea)	1	2 ea	NOS	
111	217-543	Wild Rubber Rings	3	15ea	NOS	
112	251-529	Wild MPS reticle for eyepiece	1	47	NOS	
113	<b>273-657</b>	<b>Wild Rack</b>	<b>1</b>	<b>175</b>	<b>NOS</b>	<b>f stand serial # 97545</b>
114	<b>2B260-035-1</b>	<b>Wild Pinion</b>	<b>1</b>	<b>110</b>	<b>NOS</b>	
115	<b>306-863</b>	<b>Wild Rack</b>	<b>9</b>	<b>175 ea</b>	<b>NOS</b>	<b>f stand serial #96460/964601below 94335</b>
116	315-280	Wild Adapter for lamp carrier	1	15	NOS	
117	334-840	Wild Eyepiece guard	2	10ea	NOS	rubber
118	334-860	Wild Eyepiece Guard	2	10ea	NOS	plastic
119	<b>363-447</b>	<b>Wild Pinion</b>	<b>8</b>	<b>175 ea</b>	<b>NOS</b>	<b>f stand serial #96460/964601below 94335</b>
120	367-869	Wild Ring light holder	2	30ea	NOS	
121	367-898	Wild 0.4x attachable lens	1	195	4	For M3,5,8
122	370-789	Wild Stereo 2x aux lens	1	195	4	
123	375-800	Wild Tube adapter 33mm	1	65	NOS	
124	<b>376-738</b>	<b>Wild Pinion</b>	<b>1</b>	<b>110</b>	<b>NOS</b>	
125	377-510	Wild accessory mounting screw	1	35	NOS	for mounting bino,photo tube f M3,5,8
126	377-567	Wild Stage insert	1	25	NOS	
127	383-978	Wild .05x lens	1	195	NOS	
128	389-100	Wild metal plate	1	75	NOS	
129	398-141	Wild lens	1	45	NOS	

	A	B	C	D	E	F
1	Cat Num	Description	QTY	Price	Shape	Notes
130	531-127-312	Wild Ball bearing	2	15ea	NOS	
131	706-752	Wild Ring	1	5	NOS	
132	706-753	Wild Ring	1	5	NOS	
133	810-127	Wild Knob	1	125	NOS	f stand serial #96460/964601below 94335
134		Wild Fuse holder cap	1	1	NOS	Type MTR 23 Transformer
135		Wild Pinion threads on one end	1	110	NOS	
136		Wild pinion similar to 2B260-035-1	1	110	NOS	
137		Wild Rack	1	125	NOS	
138		Wild Rack	1	125	NOS	
139		Wild Rack & pinion matched	1	175	NOS	
140		Wild Pinion	2	110 ea	NOS	
141		Wild Pinion	1	110	NOS	
142		Wild Rack	1	175	NOS	
143		Wild Slide holder screws	4	1 ea	NOS	
144		Wild nosepiece	1	175	NOS	
145		Wild metal lamphouse mount	1	25	NOS	
146		Wild .8x 3 1/4 x 4 1/4 camera cone	1	50	4	
147				11330		
148						