

	A	B	C	D	E	F
1	Cat Num	Description	QTY	Price	Shape	Notes
2		Optiphot 1 accessories				
3	78668	Nikon 100w tungsten halogen socket w/plug	1	\$75	NOS	f/87500 HMX-2/76576 PS
4		Nikon replacement lamp house lens f 12v 50w LH	1	95	NOS	Screws into stand accepts LH
5		Nikon Bulb centering target w/ND2 filter	2	20ea	NOS	Fits over field diaphragm
6		Nikon Circular stage insert	1	15	NOS	54mm dia 25mm ID
7	79084	Nikon 45mm filter holder single	4	10ea	NOS	Half moon w/handle f/ Diaphot also
8		Nikon Filter slider incl/ND32/open	2	40ea	NOS	f/ 79534 episcopic BF/DF ill Opti M
9		Nikon Optiphot instruction manual 29 pages	1		NOS	See separate listing
10	?????	Nikon Optiphot repair parts service manual			NOS	See separate listing
11		Optiphot-2 accessories				
12		Nikon plastic nosepiece cover	1	15	NOS	brand new
13		Labophot 1/2 Optiphot 1/2 Accessories				
14	79232	Nikon Cytodrive stage w/o controller	2	125	NOS	Brand new could use for stage parts
15	79443	Nikon Rubber fine focus knob	1	10	NOS	Brand new
16	84204	Nikon Tilt binocular head for XY stands	1	275	NOS	uses 79040/84220 eyepieces
17	84219	Nikon Ultra wide Binohead EP reduction ring	7	25pr	NOS	reduces diameter to use CF/CFN EP
18	84237	Nikon White condenser cap f/abbe cond	5	10ea	NOS	Brand new
19	85999	Nikon Microscope case plastic (slight use)	2	150ea	NOS	Keys included holds all access
20		Nikon Slide holder f/79010 stage (brand new)	3	45ea	NOS	Tandem geared short finger
21		Nikon Stage holding screw f/79010	1	15	NOS	Brand new unused
22		Nikon Labophot instruction manual	1		NOS	See separate listing
23		Microphot/SA Accessories				
24	79443	Nikon Rubber fine focus knob clear	5	20pr	NOS	f/SA/FXA 30mm ID approx
25	84050	Nikon Epi-U2 BF/DF/Fluorescence Illuminator	1	375	NOS	In orig box f/78648/78652 type cubes
26	84129	Nikon .8x adapter f/rear port	1	90	NOS	f/Microphot SA (new in box)
27	84132	Nikon 2.0x lens	1	90	NOS	f/Microphot SA (new in box)
28	87510	Nikon Motorstage w/o controller	1	195	NOS	Brand new could use for parts
29	87511	Nikon scanning substage f/above	1	75	NOS	Brand new incl/cond holder
30		Nikon Microphot/SA/FXA instruction manual	1		NOS	See separate listing
31						
32		Nomarski DIC for Microphot/FX/Optiphot 1				
33		<i>Contact me directly for greater detail</i>				

	A	B	C	D	E	F
1	Cat Num	Description	QTY	Price	Shape	Notes
34	84065	Nikon Nomarski DIC A cond 1.40 NA only	1	900	NOS	f/ Optiphot 1/Microphot/FX
35		DIC 20x,40x,100x obj. prisms PH1/PH3 annuli				
36	84090	Nikon Nomarski DIC prism slider only	1	800	NOS	Missing T intermediate tube
37		Missing Intermediate tube to make complete				
38		Need 79188 T intermediate tube f/ Optiphot 1				
39		Need 84016 T intermediate tube f/ Microphot				
40		Don't need 84016 if inserted into Microphot				f/Epi fluorescence attachment
41		Other DIC accessories				
42	78944	Nikon PH3 phase annulus	1	65	NOS	f/78993 TMD condenser
43	84041	Nikon DIC 10x prism f/DIC condenser	1	400	NOS	
44	84096	Nikon PH4 phase annulus	1	75	NOS	f/84047/84065 condenser
45	?????	Nikon 1/20 wave plate w/pre-polarizer graduated	1	250	NOS	f/ Optiphot 2
46						
47		Labophot 1/2 Optiphot 1/2 Fluorescence set				
48	78673	Nikon DS-Epi-FI-2 XY 2 fluor illuminator	1		NOS	
49	78648	Nikon B-2A fluor filter cube	1		NOS	
50	78652	Nikon G-2A fluor filter cube	1		NOS	
51	87530	Nikon HMX3 lamp house	1		NOS	
52	87505	Nikon Epi Quartz collector lens	1		NOS	f/ HMX2 or HMX-3
53	78692	Nikon 75w XBO Xenon socket	1		NOS	Includes free 75w Xenon merc bulb
54	78690	Nikon 75w XBO Xenon Pwr sup	1		NOS	
55	78682	Osram 75w XBO Xenon bulb	1			
56		Set price Complete		\$800		
57		Labophot 1/2 Optiphot 1/2 Fluorescence access				
58	78640	Nikon UV-1A fluor cube	2	150ea	NOS	In great shape no pinholes
59	78587	Nikon HMX Lamphouse w/collector lens	2	75ea	NOS	f/parts
60	78588	Nikon HMX HBO 100w lampsocket	1	35	NOS	f/parts broken ceramic bulb holder
61	87500	Nikon HMX-2 Lamphouse (fits below)	1	150	NOS	f/HBO/QH/XBO (good shape)
62	78592	Nikon HMX-2 HBO 100w Lamp socket fits above	1	125	NOS	2 plug igniter type no PS
63	78655	Nikon XY Epi Fluorescence Illuminator	1	150	NOS	Darker cream color
64	86966	Nikon Pinhole turret diaphragm cassette	1	295	NOS	.1/.2/.5/1/2/5mm pin holes/open/shut
65	90022	Nikon Fiber optic adapter w/f mount/col lens	1	275	NOS	Centerable fiber adjust very nice

	A	B	C	D	E	F
1	Cat Num	Description	QTY	Price	Shape	Notes
66	90021	Nikon Liquid fiber optic single f/above	1	200	NOS	12mm OD LH/ 8mm to 90022
67		Nikon auxiliary barrier slider f/ 78673	5	10ea	NOS	Plastic brand new
68		Nikon QuadFluor plastic dummy cube	3	5ea	NOS	Protects eyes from UV
69		Nikon QuadFluor right side cover	1	20	NOS	Metal keeps dust/stray light out
70		Nikon Epi Fluor light trap	1	10	NOS	Fits in place of cond cuts ref off glass
71		Nikon Fluorescence lamp centering tool	6*	35ea	NOS	Screws into nosepiece 160mm size
72		Nikon Fluorescence lamp centering tool	2*	35ea	NOS	Screws into nosepiece 160mm size
73		Diaphot/TMD Accessories				
74		Nikon Nose piece adapter	1	95	NOS	Adapter accepts fiber bundle
75		Nikon Stage insert 118mm dia	2	80ea	NOS	20mm opening brand new
76		Nikon Stage insert 118mm dia	2	80ea	NOS	10mm opening brand new
77		Nikon Stage insert 118mm dia	1	80	NOS	40mm opening brand new
78	78991	Nikon Motorized stage w/o controller	1	125	NOS	Brand new could use stage top
79	ME2	Nikon PHI phase ring insert f/78924 ELWD cond	1	50	NOS	f/4x phase objective
80		Nikon 45mm filter slider	3	25ea	NOS	120mm X 50mm Sq ea slider holds 2
81		Nikon Diaphot instruction manual	1		NOS	See separate listing
82		Nikon Diaphot repair parts service manual			NOS	See separate listing
83		Nikon Objectives 160mm				
84	78771	Nikon CF M Fluor 40x oil POL 210/0	1	\$395	NOS	New never used
85	79076	Nikon 40x Fluor .85 NA 160/.11-.23 CC	1	375	NOS	Nice clean image
86	85005	Nikon 100x Fluor 1.3 oil 160/,17 w/iris	1	425	NOS	very nice little use
87	85052	Nikon 20x Plan Ach .4 NA 160/1.2 ELWD PH2 BM	1	310	NOS	Best f/ Diaphot in original package
88	85053	Nikon 40x N Plan .55 NA 160/0-2.5 ELWD PH3 DL	1	525	NOS	Best f/ Diaphot very little use
89	85499	Nikon 100x 1.25 NA 160/.17 PH4 DL	1	370	NOS	Works fine
90	85502	Nikon 20x E DL Plan .4 NA 160/1.2 PH 2 DL	1	95	NOS	Works fine
91		Nikon Eyepieces				
92	79039	Nikon CFW 8x Eyepiece w/reticle holder	2	60 ea	NOS	Little use
93	79040	Nikon CFW 10x eyepiece w/reticle holder	1	50	NOS	Little use w/rubber eyeguard
94	79041	Nikon CFW 10x CR eyepiece w/crosshair	1	100	NOS	Focusable w/concentric circles
95	79045	Nikon CFW 10x M eyepiece w/photo reticle	1	50	NOS	Brand new in box never used
96	84220	Nikon CFWN 10x/20 eyepiece w/reticle holder	2	95 ea	NOS	New style never used
97	84229	Nikon Eyeguard for CFWN 10x eyepiece	2	10ea	NOS	NOS in original box

	A	B	C	D	E	F
1	Cat Num	Description	QTY	Price	Shape	Notes
98		FX camera Accessories				
99		Nikon SBM-1 Speed flash w/o power supply	1	125	NOS	f/HFX 2 instructions incl
100		Use below item to add Nikon DSLR to trinoc				
101	79403	Nikon PFX main body manual shutter	2		NOS	Use as a DSLR adapter
102	77129	Nikon PFX shutter camera release cable	2		NOS	2360 F to C mount adapter
103	79407	Nikkon FX35 ocular viewfinder	2		NOS	Used for focusing and composing
104	79404	Nikon FX35A adapter F mount for DSLR	2		NOS	f/ modern DSLR camera
105		Need proj lens Set price each		150		
106	79410	Nikon PL 2.5x projection lens	1	200	NOS	
107	79412	Nikon PL 5x projection lens	1	200	NOS	f/ AFXII/HFXII/UFX II
108	78834	Nikon Low pwr focusing magnifier	1	150	NOS	f/ AFXII/HFXII/UFX II
109		Nikon Remote push button switch	4	20ea	NOS	f/ AFXII/HFXII/UFX II
110	79404	Nikon FX35A adapter F mount	2	95	NOS	Nikon F mount to AFXII/HFXII/UFX II
111	79407	Nikon FX35 Ocular viewfinder	1	50	NOS	f/ AFXII/HFXII/UFX II
112		Nikon Chrome ring w/threads 43mm dia 34mm	1	25	NOS	Screws into AFXII/HFXII/UFX II
113		Nikon Chrome ring w/threads 43mm dia 17mm	1	25	NOS	Screws into AFXII/HFXII/UFX II
114	79408	Nikon FX35mm manual camera back	1	95	NOS	f/ AFXII/HFXII/UFX II
115		Nikon Shutter to control box cable	2	45ea	NOS	f/ AFXII/HFXII/UFX II
116		Labophot Optiphot Repair parts				
117	25412-004	Nikon Conecting cord for multiview	1	55	NOS	Pointer cord to PS MBF12110-1012
118		Nikon nosepiece 2-B001-028-2 f/Otiphot 1	1	45	NOS	Nosepiece holder f/dovetail
119	2G551-009	Nikon Mirror mma23000-7	2	55 ea	NOS	
120	111MBA6	Nikon Coarse Shaft & Gear	2	110 ea	NOS	f/Otiphot 1 X/XP
121	111MB#110	Nikon Coarse focus shaft	6	125 ea	NOS	
122		Nikon Lens f/bottom of 79004/05/007	4	35ea	NOS	
123		Nikon Power cord extension square 4 pin type	1	15	NOS	78662 lamp socket to 79576 trans
124	229-00-345	Nikon M plug replacement for above	3	10ea	NOS	No cord plug only
125	229-00-346	Nikon F plug replacement for above	3	10ea	NOS	No cord plug only
126	79004	Nikon Xy binocular	1	75	NOS	R side prism broken rest okay parts
127		Nikon newer accessories				
128	C-PS160	Nikon Plain stand only (focus mount below)	1	145	NOS	Post size 32mm incl b-w stage plate
129	C-FMC	Nikon Focus mount coarse and fine focus	1	175	NOS	f/ any 76mm stereomicroscope pod

	A	B	C	D	E	F
1	Cat Num	Description	QTY	Price	Shape	Notes
130	71712	Nikon MM video head Industrial	1	175	NOS	Accepts standard thread Nik objective
131	76178	Nikon Stage clips	10	15pr	NOS	In original un opened Nikon package
132		Nikon SKE/LKE/Apophot/SKT black accessories				
133	76040	Nikon Polarizing set f/ transmitted	1	175	NOS	P 89mm OD mnt a 47mm ID slip fit
134	76951	Nikon Filar eyepiece micrometer 10x scale	1	80	NOS	10mm/.01 fits any 23mm eyetube
135	77041	Nikon Phase centering telescope standard	1	50	NOS	
136	77048	Nikon Phase centering telescope short	1	45	NOS	f/ use with teaching head
137	77140	Nikon S&M Photometer for PFM w/case	1	15	NOS	Early day light meter
138	77265	Nikon short barrel Ach 10x .3 NA RECTIFIED	1	195	NOS	Rotating wave plate
139	77785	Nikon J Att mech stage 50 x 30mm movement	1	125	NOS	f/ C or P stages S-KE L-KE
140	77790	Nikon short barrel 4x .10 NA WD 9.50mm	1	100	NOS	
141	77795	Nikon short barrel 10x .25 NA 7.1mm WD	1	100	NOS	
142		Nikon 40x Phase ring insert f/M/MS inverted	1	25		
143		Phase condenser set for SKE/LKE/SKT/Apophot				
144	77040	Nikon Phase cond set incl PH tele, fitted case	1		NOS	PH rings f/ 10,20,40,100x,
145	77025	Nikon 40x DLL PH .65 NA WD .54mm	1		NOS	Short barrel objective
146	77032	Nikon 10x BM PH .30 NA 6.40 WD	1		NOS	Short barrel objective
147		Set Price		275		
148		Nikon phase cond exploded view manual			NOS	See separate listing
149		Nikon L-Ke instruction manual 24 pgs	1	20		circa 1972 PDF scan only available